






The Liberty Healthcare TILT range is specifically designed for high care patients with carers in mind. Featuring Hospital-Grade, European engineered electrics and our market leading quality steel frame.

The TILT adjustable bed features TRENDELENBURG functionality which enables forward and reverse tilt and one of the highest lift capabilities on the market. Making it beneficial for any user needing to transfer at certain heights.

#### **Functions**



- Head Raise: Improves postural support and can ease symptoms of reflux, respiratory conditions, snoring and sleep apnoea.
- Leg Raise: Assists with circulation and helps relieve tired, swollen or painful limbs.
- Zero Gravity: Helps relieve pressure on the joints, reduces compression and relaxes the muscles.
- Hi-Lo: Allows the raising of the mattress surface to an ideal height.
- Tilt: Trendelenburg (Forward and Reverse Tilt).

#### **Features**

- 5 Year Electrical Warranty\*
- 25 Year Steel Frame Warranty\*\*
- Direct Drive & Twin Drive Motor's
- TRENDELENBURG Functionality (Forward and Reverse Tilt)
- Zero Gravity Functionality
- Corded Handset
- Hi-Lo Function
- Companion Bed Capabilities

#### Standard Inclusions

- 75mm Grey Locking Castors (x4)
- 50mm Bed Base Sleeping Platform
- Foot Bar (Helps prevent the mattress from slipping whilst in motion)
- Bed Base and Foot Bar Cover

### Customisation Available

- Bed Length
- Bed Width
- Bed Height



<sup>\*5</sup> year Electrical warranty \*\*25 year Steel Frame warranty



### Optional Accessories

- Assistive Technologies e.g. Sip n Puff, Jellybean Switch and Scanning Bed Controller
- Bed Switch Mount
- Gooseneck Mount
- Battery Backup
- Side Grip Handles
- Mattress Support Bars
- Bed Rails
- Bed Rail Covers

### Cosmetic Accessories

- Upholstered Head Boards
- Upholstered Foot Boards

Bed Sizes (Please Note: Single and Double beds are longer than the standard bed length)



760mm



915mm

1070mm

1370mm









### **Optional Accessories**

- Sip n Puff (Ideal for people who have limited or no motor capability to operate switch activated devices)
- Jellybean Switch (Is designed to be used in situations where a traditional wired switch is unsuitable or inconvenient)
- Scanning Bed Controller (An essential piece of equipment used with either the Sip n Puff or the Jellybean Switch; allows people with minimal movement full control over an electric bed)
- Bed Switch Mount (Sip n Puff) (Attached under the mattress it holds the switch in an accessible position)
- Gooseneck Mount with Universal Clamp (Sip n Puff)
- Battery Backup
- Side Grip Handle
- Mattress Support Bars
- Bed Rails
- Bed Rail Covers

Please Note: Side Grip Handles, Bed Rails and Mattress Support Bars cannot be used in conjunction with Upholstered Bed Frame Side Panels.



Sip n Puff Switch



Jellybean Switch



Scanning Bed Controller



Bed Switch Mount Gooseneck Mount





Side Grip Handle



Battery Backup



Mattress Support Bar (Not suitable for Narrow Single)



Bed Rail



Bed Rail Cover



<sup>\*5</sup> year Electrical warranty \*\*25 year Steel Frame warranty



Specifications	Measurement
Under Bed Clearance (Hoist Access)	175mm (Additional castor and feet options are available if greater access height is required)
Floor - Top of Bed Base Sleeping Platform (Lowest Position*)	397mm
Floor - Top of Bed Base Sleeping Platform (Highest Position*)	880mm
Height Range	483mm
Safe Working Load	Narrow Single Long   Standard Single Long   King Single - 280kg Double   Queen - 320kg
Upper Body Profile	65°
Lower Body Profile	40°
Trendelenburg/ Reverse Trendelenburg	14° Max / 13° Max



 $<sup>^{\</sup>ast}\mbox{With}$  standard 75mm Castors and standard 50mm Bed Base Sleeping Platform.